




Stainless Steel Stone Extractors



Stainless Steel Stone Extractors


The 2.8 and 3.2 French Captura and Atlas Wire Stone Extractors feature a reinforced tip on 2.5 and 3.0 French sheaths. This configuration combines strength at the tip for multiple basket deployments with optimal fluid flow through the working channel of the scope.



Atlas Wire™

STONE EXTRACTOR

The Atlas Wire basket permits the stone to slide into the basket with ease. Atlas Wire Stone Extractors feature a 4-wire construction for maximum stone retention in a flat-wire design.



Global Product Number	Order Number	Sheath Tip French Size	Sheath French Size	Length cm	Basket	Tip Length cm
Occlusive						
G32734	AE-417115	1.7	1.7	115	4-wire	
G32788	AE-428115	2.8	2.5	115	4-wire	
G32824	AE-432115	3.2	3.0	115	4-wire	
G32739	AE-445065	4.5	4.5	65	4-wire	
Filiform						
G32735	AEF-417115	1.7	1.7	115	4-wire	5
G32737	AEF-428115	2.8	2.5	115	4-wire	5
G32730	AEF-432115	3.2	3.0	115	4-wire	5
G32740	AEF-445065	4.5	4.5	65	4-wire	5

